

### Customer and Order Information

Purchase Order *	Tag	v.1123
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in

### Base Chair

#### Fuze T50n 50° Manual Tilt Wheelchair

- |                                             |                                                                              |       |          |
|---------------------------------------------|------------------------------------------------------------------------------|-------|----------|
| <input checked="" type="checkbox"/> 2601000 | Fuze T50 50° Manual Tilt Wheelchair                                          | E1161 | \$ 3,150 |
| <input checked="" type="checkbox"/> 2601008 | Narrow Frame Upgrade Reduces overall chair width to seat width plus only 4". | K0108 | \$ 200   |

### Frame Type & Weight Capacity

#### Fixed Width Frame Select a Fixed Seat Width below.

- |                                  |                                                                                               |  |        |
|----------------------------------|-----------------------------------------------------------------------------------------------|--|--------|
| <input type="checkbox"/> 2602010 | 300 lbs Weight Capacity                                                                       |  |        |
| <input type="checkbox"/> 2602011 | 400 lbs Weight Capacity Not compatible with Transit Tie Down System.                          |  |        |
| <input type="checkbox"/> 2602025 | <b>Adjustable Width Frame</b> 250 lbs Weight Capacity. Select an Adjustable Seat Width below. |  | \$ 100 |

### Frame Options

- |                                  |                                                                                                                                                                                             |  |        |
|----------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--------|
| <input type="checkbox"/> 2602222 | <b>Take-Apart Frame</b> The Take-Apart frame allows attendants to take the Fuze apart in seconds without tools. <a href="https://pdgmobility.com/take-apart">pdgmobility.com/take-apart</a> |  | \$ 600 |
|----------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--------|
- This feature reduces the weight of the heaviest piece of the wheelchair to 29 lbs or less. Perfect for transporting it in the trunk of a car or SUV.  
 Not compatible with Power Tilt, Ventilator Tray, Ventilator Battery Tray, DRB, User Wheel Locks.  
 Take-Apart rear crossbar is not width adjustable. PDG will provide one no-charge crossbar in a different width over the lifetime of each chair ordered with adjustable width frame (can be orderd with the chair or later). Note: Crossbar for the T50n is 4" shorter than the Seat Width.



- |                                  |                                       |  |  |
|----------------------------------|---------------------------------------|--|--|
| <input type="checkbox"/> 2602300 | <b>Transit Tie Down System (TTDS)</b> |  |  |
|----------------------------------|---------------------------------------|--|--|
- TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.  
 Tested to max. Weight Capacity of 300 lbs.

### Seat Width

Overall chair width is +4.0" for 12" RW

#### Fixed Seat Widths Order with Fixed Frame.

- |                                  |                |                                  |                |       |        |
|----------------------------------|----------------|----------------------------------|----------------|-------|--------|
| <input type="checkbox"/> 2603016 | 16" Seat Width | <input type="checkbox"/> 2603017 | 17" Seat Width |       |        |
| <input type="checkbox"/> 2603018 | 18" Seat Width | <input type="checkbox"/> 2603019 | 19" Seat Width |       |        |
| <input type="checkbox"/> 2603020 | 20" Seat Width | <input type="checkbox"/> 2603021 | 21" Seat Width | E2201 | \$ 375 |
| <input type="checkbox"/> 2603022 | 22" Seat Width | <input type="checkbox"/> 2603023 | 23" Seat Width | E2201 | \$ 375 |
| <input type="checkbox"/> 2603024 | 24" Seat Width |                                  |                | E2202 | \$ 575 |

#### Adjustable Seat Widths Order with Adjustable Frame.

- |                                  |                            |                                  |                            |       |        |
|----------------------------------|----------------------------|----------------------------------|----------------------------|-------|--------|
| <input type="checkbox"/> 2603216 | 16"-19" Range - Set at 16" | <input type="checkbox"/> 2603217 | 16"-19" Range - Set at 17" |       |        |
| <input type="checkbox"/> 2603218 | 16"-19" Range - Set at 18" | <input type="checkbox"/> 2603219 | 16"-19" Range - Set at 19" |       |        |
| <input type="checkbox"/> 2603319 | 19"-20" Range - Set at 19" |                                  |                            |       |        |
| <input type="checkbox"/> 2603320 | 19"-20" Range - Set at 20" |                                  |                            | E2201 | \$ 375 |



# FUZE T50n

## Order Form - USA

Prices in U.S. Dollars. January 1, 2024

### Seat Depth

#### Standard Seat Depths (Adjustable)

<input type="checkbox"/> 2604015	15" Seat Depth	<input type="checkbox"/> 2604016	16" Seat Depth	<input type="checkbox"/> 2604017	17" Seat Depth		
<input type="checkbox"/> 2604018	18" Seat Depth	<input type="checkbox"/> 2604019	19" Seat Depth				
<input type="checkbox"/> 2604020	20" Seat Depth					E2203	\$ 500

#### Custom Seat Depths

Requires Stroller Handle Back - Straight. Not compatible with Reclining or Fold-down Back.

<input type="checkbox"/> 2604021	21" Seat Depth			E2203	\$ 525
<input type="checkbox"/> 2604022	22" Seat Depth			E2204	\$ 800



### Frame Color

#### PDG Premium Colors

View color swatches on PDG's website. [pdgmobility.com/colors](https://pdgmobility.com/colors)

<input type="checkbox"/> 2605001	Black Shimmer	<input type="checkbox"/> 2605015	Sky Blue	<input type="checkbox"/> 2605016	Poppy Red
<input type="checkbox"/> 2605003	Fine Texture Black	<input type="checkbox"/> 2605010	Lollypop Blue	<input type="checkbox"/> 2605011	Lollypop Red
<input type="checkbox"/> 2605017	Arctic Grey	<input type="checkbox"/> 2605018	Tropical Teal	<input type="checkbox"/> 2605019	Tangerine Orange
<input type="checkbox"/> 2605020	Olive Green	<input type="checkbox"/> 2605021	Lush Lilac	<input type="checkbox"/> 2605012	Bubble Gum Pink
<input type="checkbox"/> 2605022	Leaf Green	<input type="checkbox"/> 2605023	Boysenberry Purple	<input type="checkbox"/> 2605014	Daffodil Yellow
<input type="checkbox"/> 2605200	Custom Color	Contact PDG for specifications.			\$ 300

### Armrests

<input type="checkbox"/> 2606000	No Armrests				
<b>Height Adjustable Armrests (Single Post)</b>				E0973	\$ 345
<input type="checkbox"/> 2606114	10"-14" Range	<input type="checkbox"/> 2606116	12"-16" Range		
<input type="checkbox"/> 2606100	Fixed Height Armrests (12")	Full Length Armpad only.			
	<b>Armrest Height</b>	If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.			
	Seat Pan to top of Armrest:		Does not include height of Armpad		
	<b>Armrest Position</b>	Indicate preference for Armrest position. If left blank, the default setting will be used.			
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.					
<input type="checkbox"/> 2606402	Mid	<input type="checkbox"/> 2606403	Front		
<b>Locking Flip-up Armrests</b>				E0973	\$ 345
Not recommended for users with excessive strength or tone. Not compatible with Fold-down Back.					
<input type="checkbox"/> 2606310	Set at 10"	<input type="checkbox"/> 2606311	Set at 11"	<input type="checkbox"/> 2606312	Set at 12"
<input type="checkbox"/> 2606313	Set at 13"	<input type="checkbox"/> 2606314	Set at 14"		

### Armpads

<b>Regular Armpads</b>					
<input type="checkbox"/> 2607100	Desk Length - 2" x 10"	<input type="checkbox"/> 2607101	Full Length - 2" x 14"		
					\$ 115
<b>Waterfall Armpads</b>					
<input type="checkbox"/> 2607105	Desk Length - 2" x 10"	<input type="checkbox"/> 2607106	Full Length - 2" x 14"		
					\$ 200
<b>Wide Armpads</b> May interfere with wheels in some tilt configurations.					
<input type="checkbox"/> 2607110	Desk Length - 3" x 10"	<input type="checkbox"/> 2607111	Full Length - 3" x 12"		
					\$ 135
<b>Gel Armpads</b> May interfere with wheels in some tilt configurations.					
<input type="checkbox"/> 2607120	Full Length - 2" x 12"	<input type="checkbox"/> 2607121	Full Length - 3.5" x 12"		
				E2209	\$ 125
<b>Contour Arm Trough with Hand Pad</b>					
Not compatible with Lap Tray, Locking Flip-up Armrests. May interfere with wheels in some tilt configurations.					
<input type="checkbox"/> 2607130	Right (4.5" x 18")	<input type="checkbox"/> 2607131	Left (4.5" x 18")		

### Armrest Options

<b>Armrest Wideners</b> Move Armrests outward by up to 1" on each side. Seat size does not change.					
<input type="checkbox"/> 2607210	1" Total Additional Width				\$ 95
<input type="checkbox"/> 2607211	2" Total Additional Width				\$ 195



# FUZE T50n


## Order Form - USA

Prices in U.S. Dollars. January 1, 2024

### Seat Pan

<input type="checkbox"/> 2608000	<b>No Seat Pan</b>		
<input type="checkbox"/> 2608001	<b>Solid Aluminum Seat Pan</b>	E2231	\$ 165

### Seat Height

 Don't forget

Select desired Seat-to-Floor Height    STF Height is adjustable.


12" Rear Wheels					
Caster size	16"	17"	18"	19"	20"
5"	<input type="checkbox"/>	<input type="checkbox"/>			
6"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
8"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
6" Soft Roll	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

### Casters

Urethane Casters										
<input type="checkbox"/>	2611051	5" x 1.25" Caster size	<input type="checkbox"/>	2611061	6" x 1.25" Caster size	<input type="checkbox"/>	2611081	8" x 1.25" Caster size		
<input type="checkbox"/>	2611063	6" x 1.50" Caster size							\$	85
<input type="checkbox"/>	2611084	8" x 2.00" Caster size							\$	135
Pneumatic Casters										
<input type="checkbox"/>	2611161	6" x 1.25" Caster size						E2214 & E2215	\$	110
<input type="checkbox"/>	2611181	8" x 1.25" Caster size						E2214 & E2215	\$	120
<input type="checkbox"/>	2611184	8" x 2.00" Caster size						E2214 & E2215	\$	150
<input type="checkbox"/>	2611260	<b>Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters</b>							\$	275

### Rear Wheels

☒ 2612012    **12" Composite Wheels**



**Rear Wheel Position**    Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.  
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

☐ 2612401    Rear    More stable    ☐ 2612402    Mid    ☐ 2612403    Front    Less stable

### Axle Type

<input type="checkbox"/> 2613000	<b>Quick Release Rear Axles</b>
<input type="checkbox"/> 2613010	<b>Fixed Rear Axles</b>

### Rear Tires

<input type="checkbox"/> 2614000	<b>Urethane Tires</b>		
<input type="checkbox"/> 2614010	<b>Pneumatic Tires</b> May result in a finished Seat Height 1/2" greater than stated.	E2211 & E2212	\$ 60
<input type="checkbox"/> 2614020	<b>Pneumatic Tires with Flat-Free Inserts</b> May result in a finished Seat Height 1/2" greater than stated.	E2211 & E2213	\$ 95

### Handrims

<input checked="" type="checkbox"/> 2615000	<b>No Handrims</b>
---------------------------------------------	--------------------

### Wheel Locks & Brakes

<b>User Wheel Locks</b> Not compatible with Take-Apart Frame.		\$ 150
<input type="checkbox"/> 2616010    Push to Lock	<input type="checkbox"/> 2616020    Pull to Lock	
<input type="checkbox"/> 2616030 <b>Attendant Wheel Locks (Push to Lock)</b>		
<input type="checkbox"/> 2616005 <b>Foot Operated Wheel Lock</b> Rear-mounted unilateral single-activation lock. Not compatible with Ventilator Tray.		\$ 595
<input type="checkbox"/> 2616300 <b>Drum Brakes</b> Available with Urethane Tires only. Compatible with 17.5", 18.5", 19.5" STF Height only.		<b>COMING SOON</b> E2228    \$ 795
<b>User-activated Unilateral Wheel Lock</b> Pull to Lock. Available up to 22" Seat Width only.		\$ 795
<input type="checkbox"/> 2616310    Left Activation	<input type="checkbox"/> 2616311    Right Activation	



# FUZE T50n

## Order Form - USA

Prices in U.S. Dollars. January 1, 2024

### Wheel Lock Options

<input type="checkbox"/>	2616100	Wheel Lock Extensions	E0961	\$	80
--------------------------	---------	-----------------------	-------	----	----

### Anti-Tippers

<input checked="" type="checkbox"/>	2617010	Rear Anti-Tippers	Required Safety Option.	E0971	\$	125
-------------------------------------	---------	-------------------	-------------------------	-------	----	-----

### Back Posts

#### Back Post Style

<input type="checkbox"/>	2618100	Stroller Handle Back	<input type="checkbox"/>	2618105	Cane Handle Back	Not compatible with Headrest Mounting Plate.
--------------------------	---------	----------------------	--------------------------	---------	------------------	----------------------------------------------

#### Back Post Type

<input type="checkbox"/>	2618200	Straight	<input type="checkbox"/>	2618205	8° Bend
--------------------------	---------	----------	--------------------------	---------	---------

#### Back Post Height

 Measured from the Seat Pan to the Handle.

<input type="checkbox"/>	2618420	Short (17")	<input type="checkbox"/>	2618425	Tall (23")
--------------------------	---------	-------------	--------------------------	---------	------------

#### Back Post Folding

 Fold-down Back is not compatible with Dynamic Rocker Back Interface, Reclining Backrest.

<input type="checkbox"/>	2619220	Fold-down	<input type="checkbox"/>	2619221	Non-folding
--------------------------	---------	-----------	--------------------------	---------	-------------

#### Seat to Back Angle

<input type="checkbox"/>	2619090	Set at 90°	<input type="checkbox"/>	2619096	Set at 98°	<input type="checkbox"/>	2619102	Set at 106°
<input type="checkbox"/>	2619108	Set at 114°	<input type="checkbox"/>	2619114	Set at 122°			

### Back Options

#### Angle/Height Adjustable Stroller Handle

 Not compatible with Headrest Mounting Plate or Cane Handle Back.

<input type="checkbox"/>	2619010	Button on Outside	<input type="checkbox"/>	2619011	Button on Inside	\$	150
--------------------------	---------	-------------------	--------------------------	---------	------------------	----	-----

#### Removable Cane Handle Push Bar

 Not compatible with Stroller Handle Back.

<input type="checkbox"/>	2619050					\$	140
--------------------------	---------	--	--	--	--	----	-----

#### Endomorph Back

 Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12".

<input type="checkbox"/>	2618102	1" Offset per side	<input type="checkbox"/>	2618103	1.5" Offset per side	<input type="checkbox"/>	2618104	2" Offset per side
--------------------------	---------	--------------------	--------------------------	---------	----------------------	--------------------------	---------	--------------------

#### Headrest Mounting Plate

 Not compatible with Angle/Height Adjustable Stroller Handle, Cane Handle Back. Not required if mounting to a solid Back.

<input type="checkbox"/>	2619212					\$	95
--------------------------	---------	--	--	--	--	----	----

#### Seating Dynamics Dynamic Rocker Back Interface

E2398	\$	2,995
-------	----	-------

Not compatible with Take-Apart Frame, Fold-down Back, 30° Recline, Standard Upholstery. Only available with Stroller Handle - 8° Bend.  
Finished Back Post Height is 7" taller than selected Back Post Height. Maximum compatible Seat Depth is 20". 90° and 98° Seat Angle only.

### Back Upholstery

<input type="checkbox"/>	2620000	No Back Upholstery			
--------------------------	---------	--------------------	--	--	--

#### Standard Upholstery

 Not available with Cane Handle Back.

E2611 / E2612	\$	195
---------------	----	-----

<input type="checkbox"/>	2620418	Short (18" Height)	<input type="checkbox"/>	2620421	Regular (21" Height)
--------------------------	---------	--------------------	--------------------------	---------	----------------------

#### Tension Adjustable Upholstery

 Height of Upholstery is 16". Not available with Cane Handle Back.

E2611 / E2612	\$	295
---------------	----	-----

#### Position of Upholstery

 Select the desired position of the top of the Upholstery relative to the Seat Pan.

17" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>								

23" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Regular					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Tension Adj				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

### Tilt Activation

#### Attendant Tilt only

 Select trigger mounting location below.

<input type="checkbox"/>	2623011	Right Back Post	<input type="checkbox"/>	2623012	Left Back Post
--------------------------	---------	-----------------	--------------------------	---------	----------------

See next page for more options >



# FUZE T50n

## Order Form - USA

Prices in U.S. Dollars. January 1, 2024

<b>Power Tilt</b> 40° tilt range. Includes 2 batteries for off-chair charging. <a href="https://pdgmobility.com/power-tilt">pdgmobility.com/power-tilt</a>		E1002	
<input type="checkbox"/> 2623400	Standard Power Tilt (250 lbs Weight Capacity)	\$ 2,295	
<input type="checkbox"/> 2623401	Standard Power Tilt (350 lbs Weight Capacity) Requires Fixed Frame.	\$ 2,595	
<input type="checkbox"/> 2623410	Power Tilt with Manual Tilt Override (250 lbs Weight Capacity) Allows for instant manual release of tilt in case of power tilt failure.	\$ 2,795	
<input type="checkbox"/> 2623411	Power Tilt with Manual Tilt Override (350 lbs Weight Capacity) Allows for instant manual release of tilt in case of power tilt failure. Requires Fixed Frame.	\$ 3,095	
<div>Select Control Switch</div>	Toggle Switch	<input type="checkbox"/> 2623421 Installed on Right side	<input type="checkbox"/> 2623422 Installed on Left side
	Button Switch	<input type="checkbox"/> 2623431 Installed on Right side	<input type="checkbox"/> 2623432 Installed on Left side
	Stealth Egg Switch	<input type="checkbox"/> 2623440 Not installed	
	Optional for Egg Switch:	<input type="checkbox"/> 2623445 Swing Away Mounting Hardware	Not compatible with all Headrest makes and models.
Specify Headrest to be used: _____		E1028 \$ 625	
<b>Tilt Limit</b> Optional. If no tilt limit is selected, the full tilt range will be available.			
<input type="checkbox"/> 2623220	Set at 20°	<input type="checkbox"/> 2623230 Set at 30°	<input type="checkbox"/> 2623240 Set at 40°

### Recline

<input type="checkbox"/> 2623310	<b>Reclining Backrest (Manual)</b>	E1225 \$ 895
30° maximum recline. Not compatible with Cane Handle Back, Fold-down Back. Not available with 22" Seat Depth. Finished Back Post Height is 3" taller than selected Back Post Height.		

### Front Rigging

<input type="checkbox"/> 2625000	<b>No Front Rigging</b>	<input type="checkbox"/> 2625001 <b>Sector Blocks only</b>
<b>Easy-Out™ Front Rigging</b> Effortless insertion and removal. No Hanger Bracket remains on frame. <a href="https://pdgmobility.com/easy-out">pdgmobility.com/easy-out</a>		
<input type="checkbox"/> 2625410	Fixed 70°	\$ 95
<input type="checkbox"/> 2625420	Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°)
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.		
<b>Swing Away Front Rigging</b> Legacy Design. Not compatible with Easy-Out. Requires pre-tilt of 5°.		
<input type="checkbox"/> 2625130	Swing Away 60°	<input type="checkbox"/> 2625140 Swing Away 70°
<input type="checkbox"/> 2625160	Swing Away 90°	Includes 11"-16" Extension Tubes. Not compatible with 8" Casters.
<input type="checkbox"/> 2625200	Swing Away Elevating Legrests	Includes Calf Pads.
<input type="checkbox"/> 2625300	1" Offset Hangers for Swing Away	E0990 \$ 415
<input type="checkbox"/> 2625110	<b>Standard Plug In 70° Front Rigging</b>	Legacy Design. Not compatible with Easy-Out.

### Footplates

<input type="checkbox"/> 2626010	<b>Composite Footplates</b>	
<input type="checkbox"/> 2626020	<b>Aluminum Angle Adjustable Footplates</b>	K0040 \$ 200
<input type="checkbox"/> 2626030	<b>One Piece Fixed Footplate</b>	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.
<input type="checkbox"/> 2626040	<b>One Piece Flip-up Footplate</b>	Requires selection of Composite or Aluminum Angle Adjustable.
<input type="checkbox"/> 2626050	<b>Interlocking Angle Adjustable Footplates</b>	Available for Seat Widths up to 21" only.

### Footrest/Legrest Accessories

<b>Calf Pad</b> Not compatible with Gel Padding.		
<input type="checkbox"/> 2627012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.
<input type="checkbox"/> 2627010	For Swing Away Front Rigging	
<input type="checkbox"/> 2627020	<b>Gel Padding</b>	Not compatible with Calf Pad or Easy-Out Front Rigging.
<b>Gel Knee Pad</b>		
<input type="checkbox"/> 2627032	For Easy-Out™ Angle Adjustable Front Rigging	
<input type="checkbox"/> 2627030	For Swing Away ELR Front Rigging	
<input type="checkbox"/> 2627040	<b>Length Adjustable Heel Loops</b>	E0951 \$ 50
<input type="checkbox"/> 2627050	<b>Bodypoint Padded Calf Strap</b>	Sizes greater than 22" are not padded.



# FUZE T50n

## Order Form - USA

Prices in U.S. Dollars. January 1, 2024

### Front Rigging Extension Tubes

#### Regular Extension Tubes

No selection required for 90° Front Rigging

- ☐ 2628010 Short 14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging
- ☐ 2628020 Long 18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging

\$ 75



#### Lower Leg Length

If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.

Seat Pan to Foot Plates: \_\_\_\_\_ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.

### Miscellaneous Options

#### Stealth Flip Down Abductor

Requires Solid Seat Pan.

E0957 & E1028 \$ 650

- ☐ 2630060 Small Pad 3"x4"x3"x2" size ☐ 2630062 Large Pad 5"x5"x4"x3" size

#### QR Bracket - Accessory Mount

E1028 \$ 200

#### Gel Hip Guides

Reduces available Seat Width by 2".

E0956 & E1028 \$ 495

#### Hip Belt Mount

Provides improved mounting options for Bodypoint Hip Belt or other third-party belts.

E1028 \$ 55

#### Ventilator Tray

Reduces ground clearance. Chair may be unstable when unoccupied. Minimum Seat Width is 20".  
Not compatible with Power Tilt, Foot Operated Wheel Locks.

E1029 \$ 295

#### Ventilator Bracket

Compatible with LTV1 (laptop) ventilator.

E1029 \$ 295

#### Clear Lap Tray (Slide On)

Available for Seat Widths of 17"-22" only. Not compatible with Contour Arm Trough.

E0950 \$ 490

#### Q'Straint Lap Belt

Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length.

\$ 245

#### Wheelchair Backpack

Black with PDG Logo. 17"x14" size.

\$ 145

#### Adhesive Back Hook Strip, 2" x 30"

For Cushion mounting.

\$ 10

### On-Chair Seating & Positioning

Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. [pdgmobility.com/seating](https://pdgmobility.com/seating)

### Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions. Use of HCPCS codes does not guarantee coverage or payment.

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.