



Customer and Order Information

Purchase Order *	Tag	v.0901
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs Height _____ ft/in	* required

Base Chair

<input checked="" type="checkbox"/> 1201000	Stellar 45° Manual Tilt Wheelchair	E1161	\$ 3,150
<input type="checkbox"/>	Also available: the lightweight Stellar GLT with aluminum frame. See separate order form. pdgmobility.com/glt		
<input type="checkbox"/>	Choose the Stellar Impact for highly agitated users and the Stellar HD for bariatric users. See separate order forms.		

Weight Capacity

Chair Weight Capacity			
<input type="checkbox"/> 1202000	250 lbs Capacity	Available for 14" to 22" Seat Width.	
<input type="checkbox"/> 1202050	300 lbs Capacity	Available for 14" to 22" Seat Width.	\$ 200
<input type="checkbox"/> 1205006	450 lbs Capacity	Available for 16" to 24" Seat Width. Requires 20" Composite or 24" Heavy Duty Spoke Wheels. Fixed Axles only. Not compatible with TTDS.	\$ 500

Frame Options

<input type="checkbox"/> 1202200	Extended Base with Rear Floor Jacks	K0108	\$ 900
Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant Wheel Lock, HD Anti-Tipper, IV Pole, O2 Holder.			
<input type="checkbox"/> 1202201	Shortened Strut Bracket	Increases heel clearance for foot propelling. Available only for 23" and 24" Seat Width.	\$ 125
<input type="checkbox"/> 1203300	Transit Tie Down System (TTDS)	TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible Seat Widths greater than 22". Tested to max. Weight Capacity of 300 lbs.	

Seat Width

Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

Standard Seat Widths			
<input type="checkbox"/> 1203014	14" Seat Width	<input type="checkbox"/> 1203015	15" Seat Width
<input type="checkbox"/> 1203017	17" Seat Width	<input type="checkbox"/> 1203018	18" Seat Width
<input type="checkbox"/> 1203020	20" Seat Width	<input type="checkbox"/> 1203021	21" Seat Width
<input type="checkbox"/> 1203016	16" Seat Width	<input type="checkbox"/> 1203019	19" Seat Width
<input type="checkbox"/> 1203022	22" Seat Width	E2201	\$ 375
Custom Seat Widths Requires 450 lbs Weight Capacity Frame. Extended delivery lead times may apply.			
<input type="checkbox"/> 1203023	23" Seat Width	E2201	\$ 375
<input type="checkbox"/> 1203024	24" Seat Width	E2202	\$ 575

Seat Depth

Standard Seat Depths (Adjustable)			
<input type="checkbox"/> 1204016	16" Seat Depth	<input type="checkbox"/> 1204017	17" Seat Depth
<input type="checkbox"/> 1204018	18" Seat Depth	<input type="checkbox"/> 1204019	19" Seat Depth
<input type="checkbox"/> 1204020	20" Seat Depth	E2203	\$ 500
Custom Seat Depths			
<input type="checkbox"/> 1204021	21" Seat Depth	E2203	\$ 525
<input type="checkbox"/> 1204022	22" Seat Depth	E2204	\$ 800



Frame Color

PDG Premium Colors

- | | | |
|--|--|---|
| <input type="checkbox"/> 1205000 Tracer Red | <input type="checkbox"/> 1205001 Gloss Black Shimmer | <input type="checkbox"/> 1205002 Silver |
| <input type="checkbox"/> 1205003 Fine Texture Black | <input type="checkbox"/> 1205004 Orange Fusion | <input type="checkbox"/> 1205005 Karl's Imagination |
| <input type="checkbox"/> 1205006 Blue Shimmer | <input type="checkbox"/> 1205007 Gloss White | <input type="checkbox"/> 1205008 Dark Titanium |
| <input type="checkbox"/> 1205009 Geyser Green | <input type="checkbox"/> 1205010 Lollipop Blue | <input type="checkbox"/> 1205011 Lollipop Red |
| <input type="checkbox"/> 1205012 Bubble Gum Pink | <input type="checkbox"/> 1205013 Berry Purple | <input type="checkbox"/> 1205014 Daffodil Yellow |
| <input type="checkbox"/> 1205200 Custom Color Contact PDG for specifications. | | |

\$ 300

Armrests

- | | | | |
|---|--|--|--------|
| <input type="checkbox"/> 1206000 No Armrests | | | |
| Single Post Height Adjustable Armrests Not compatible with the 450 lbs Weight Capacity Frame. | | E0973 | \$ 295 |
| <input type="checkbox"/> 1206111 8"-11" Range | <input type="checkbox"/> 1206114 10"-14" Range | | |
| Dual Post Height Adjustable Armrests Required for the 450 lbs Weight Capacity Frame. | | E0973 | \$ 295 |
| <input type="checkbox"/> 1206211 8"-11" Range | <input type="checkbox"/> 1206214 10"-14" Range | <input type="checkbox"/> 1206216 12"-16" Range | |
| Fixed Height Armrests (12") Full Length Armpad only. | | | |
| <input type="checkbox"/> 1206100 Single Post | <input type="checkbox"/> 1206200 Dual Post | | |
|  Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.
Seat Pan to top of Armrest: _____ Does not include height of Armpad | | | |
|  Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used.
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference. | | | |
| <input type="checkbox"/> 1206401 Rear | <input type="checkbox"/> 1206402 Mid | <input type="checkbox"/> 1206403 Front | |
| Locking Flip-up Armrests Not recommended for users with excessive strength or tone. | | E0973 | \$ 295 |
| Not compatible with Fold Down Back, Clothing Guards, Back Upholstery. Available with Full Length and Desk Length Armpads. | | | |
| <input type="checkbox"/> 1206310 Set at 10" | <input type="checkbox"/> 1206311 Set at 11" | <input type="checkbox"/> 1206312 Set at 12" | |
| <input type="checkbox"/> 1206313 Set at 13" | <input type="checkbox"/> 1206314 Set at 14" | | |

Armpads

- | | | | |
|--|---|-------|--------|
| Regular Armpads | | | |
| <input type="checkbox"/> 1207100 Desk Length - 2" x 10" | <input type="checkbox"/> 1207101 Full Length - 2" x 14" | | |
| Wide Armpads May interfere with wheels in some tilt configurations. | | | \$ 200 |
| <input type="checkbox"/> 1207110 Desk Length - 3" x 10" | <input type="checkbox"/> 1207111 Full Length - 3" x 12" | | |
| Gel Armpads May interfere with wheels in some tilt configurations. | | | \$ 135 |
| <input type="checkbox"/> 1207120 Full Length - 2" x 12" | <input type="checkbox"/> 1207121 Full Length - 3.5" x 12" | | |
| Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations. | | E2209 | \$ 125 |
| <input type="checkbox"/> 1207130 Right (4.5" x 18") | <input type="checkbox"/> 1207131 Left (4.5" x 18") | | |


Armrest Options

- | | | | |
|---|--|--|--------|
| Armrest Wideners Move Armrests outward by 1" on each side. Do not increase seat size. | | | |
| <input type="checkbox"/> 1207210 For Single Post Armrests May cause tilt interference at lower Seat Height and/or Seat Depth settings. | | | \$ 195 |
| <input type="checkbox"/> 1207220 For Dual Post Armrests Limits tilt. Not compatible with 20", 22", or 24" Rear Wheels, Clothing Guards. | | | \$ 375 |
| Clothing Guards Not compatible with Flip-up Armrests, Armrest Wideners. | | | |
| <input type="checkbox"/> 1207310 For Single Post Armrests Available for 10"-14" Height range only. | | | \$ 110 |
| <input type="checkbox"/> 1207320 For Dual Post Armrests | | | \$ 165 |
| Armrest Extensions Not compatible with Flip-up or Single Post Armrests. | | | \$ 75 |
| <input type="checkbox"/> 1207400 3" Extension - Forward | <input type="checkbox"/> 1207401 3" Extension - Rearward | | |

Seat Pan

<input type="checkbox"/> 1208000	No Seat Pan		
<input type="checkbox"/> 1208001	Solid Aluminum Seat Pan	E2231	\$ 135

Seat Height

 **Don't forget** **Select desired Seat-to-Floor Height** STF Height is adjustable.

6" Casters							
RW	14"	15"	16"	17"	18"	19"	20"
12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

8" Casters							
RW	14"	15"	16"	17"	18"	19"	20"
16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Casters

Urethane Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/> 1211061	6" x 1.25" Caster size	<input type="checkbox"/> 1211081	8" x 1.25" Caster size		
<input type="checkbox"/> 1211063	6" x 1.50" Caster size				\$ 85
<input type="checkbox"/> 1211084	8" x 2.00" Caster size				\$ 135

Pneumatic Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/> 1211161	6" x 1.25" Caster size	E2214 & E2215	\$ 110
<input type="checkbox"/> 1211181	8" x 1.25" Caster size	E2214 & E2215	\$ 120
<input type="checkbox"/> 1211184	8" x 2.00" Caster size	E2214 & E2215	\$ 150

☐ 1211260 **Frog Legs 6" x 1.40" Soft Roll Casters** \$ 95

☐ 1211400 **Frog Legs Suspension Forks** Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster. E1015 \$ 395

Rear Wheels


Composite Wheels

<input type="checkbox"/> 1212012	12" Wheel size	<input type="checkbox"/> 1212016	16" Wheel size	<input type="checkbox"/> 1212020	20" Wheel size
<input type="checkbox"/> 1212022	22" Wheel size	<input type="checkbox"/> 1212024	24" Wheel size		

Spoke Wheels \$ 75

<input type="checkbox"/> 1212120	20" Wheel size	<input type="checkbox"/> 1212122	22" Wheel size	<input type="checkbox"/> 1212124	24" Wheel size
----------------------------------	----------------	----------------------------------	----------------	----------------------------------	----------------

☐ 1212224 **24" Heavy Duty Spoke Wheels** \$ 350

 **Rear Wheel Position** Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

<input type="checkbox"/> 1212401	Rear	<input type="checkbox"/> 1212402	Mid	<input type="checkbox"/> 1212403	Front
----------------------------------	------	----------------------------------	-----	----------------------------------	-------

Axle Type

<input type="checkbox"/> 1213000	Quick Release Rear Axles
<input type="checkbox"/> 1213010	Fixed Rear Axles

Rear Tires

<input type="checkbox"/> 1214000	Urethane Tires		
<input type="checkbox"/> 1214010	Pneumatic Tires	May result in a finished Seat Height 0.5" greater than stated.	E2211 & E2212 \$ 60
<input type="checkbox"/> 1214020	Pneumatic Tires with Flat-Free Inserts	May result in a finished Seat Height 0.5" greater than stated. Not available with Spoke Wheels.	E2211 & E2213 \$ 95

Handrims

Handrims are available only with 20", 22", and 24" Rear Wheels.

<input type="checkbox"/> 1215000	No Handrims	
<input type="checkbox"/> 1215010	Standard Aluminum Handrims	
<input type="checkbox"/> 1215020	Plastic Coated Handrims	\$ 150

Wheel Locks & Brakes

User Wheel Locks

<input type="checkbox"/> 1216010	Push to Lock	<input type="checkbox"/> 1216020	Pull to Lock		
<input type="checkbox"/> 1216030	Attendant Wheel Locks (Push to Lock) Not compatible with TTDS or the 300 & 450 lbs Weight Capacity Frames.			\$	85
<input type="checkbox"/> 1216040	User & Attendant Wheel Locks Available with 20", 22", or 24" Rear Wheels only. Not compatible with Drum Brakes.			\$	125
<input type="checkbox"/> 1216300	Drum Brakes	Dual trigger activation (right/left side activated independently). Available with 20" Rear Wheels only.		E2228	\$ 795

Wheel Lock Options

<input type="checkbox"/> 1216100	Wheel Lock Extensions	Not compatible with Attendant Wheel Lock.	E0961	\$	60
<input type="checkbox"/> 1216200	Hill Holders	Requires Pneumatic Rear Tires with Flat-Free Inserts.	E0974	\$	450

Anti-Tippers

<input checked="" type="checkbox"/> 1217010	Rear Anti-Tippers	Required Safety Option.	E0971	\$	75
<input type="checkbox"/> 1217020	Front Anti-Tippers		E0971	\$	795
<input type="checkbox"/> 1217030	Heavy Duty Rear Anti-Tipper	Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.	E0971	\$	895

Back Posts

Back Post Style

<input type="checkbox"/> 1218010	Push Bar Back	<input type="checkbox"/> 1218020	Cane Handle Back	Seat Widths of 14"-22" only. Not compatible with Headrest Mounting Plate.	
----------------------------------	---------------	----------------------------------	------------------	---	--

Back Post Height

Measured from the Seat Pan to the Handle.

<input type="checkbox"/> 1218425	Short (22")	<input type="checkbox"/> 1218431	Tall (25")
----------------------------------	-------------	----------------------------------	------------

Seat to Back Angle

<input type="checkbox"/> 1219090	Set at 90°	<input type="checkbox"/> 1219098	Set at 98°	<input type="checkbox"/> 1219106	Set at 106°
<input type="checkbox"/> 1219114	Set at 114°	<input type="checkbox"/> 1219122	Set at 122°		

Back Options

	Angle/Height Adjustable Push Bar Handle	Not compatible with Headrest Mounting Plate or Cane Handle Back.		\$	150
<input type="checkbox"/> 1219010	Button on Outside	<input type="checkbox"/> 1219011	Button on Inside		
<input type="checkbox"/> 1219210	10" Backrest Extension	Requires the Angle/Height Adjustable Push Bar Handle.	E0966	\$	150
<input type="checkbox"/> 1219212	Headrest Mounting Plate	Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.		\$	95
<input type="checkbox"/> 1219220	Fold Down Back	Not compatible with Reclining Back.			
<input type="checkbox"/> 1219230	Spring Back Suspension	Not compatible with Reclining Back.		\$	105
<input type="checkbox"/> 1219240	Seating Dynamics Dynamic Rocker Back Interface	Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Requires Adjustable Tension Upholstery or rigid Backrest. Finished Back Post Height 3" taller than selected Back Post Height.	K0108	\$	2,995

Back Upholstery

<input type="checkbox"/> 1220000	No Back Upholstery				
	Standard Upholstery	Not available with Cane Handle Back.	E2611 / E2612	\$	195
<input type="checkbox"/> 1220418	Short (18" Height)	<input type="checkbox"/> 1220421	Regular (21" Height)	<input type="checkbox"/> 1220424	Tall (24" Height)
<input type="checkbox"/> 1220450	Tension Adjustable Upholstery	Height of Upholstery is 16". Not available with Cane Handle Back.	E2611 / E2612	\$	295

Position of Upholstery

Select the desired position of the top of the Upholstery relative to the Seat Pan.

22" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

25" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tall							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Tilt Activation

Attendant Tilt only Select trigger mounting location below.		
<input type="checkbox"/> 1223011 Right Back Post	<input type="checkbox"/> 1223012 Left Back Post	
User Tilt only Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting location below.		
<input type="checkbox"/> 1223021 Right Armrest	<input type="checkbox"/> 1223022 Left Armrest	\$ 100
Attendant & User Tilt Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting locations below.		
<input type="checkbox"/> 1223031 Right Back Post	<input type="checkbox"/> 1223032 Left Back Post	
<input type="checkbox"/> 1223041 Right Armrest	<input type="checkbox"/> 1223042 Left Armrest	
<input type="checkbox"/> 1223110 Trigger Extension		\$ 50
Tilt Limit Optional. If no tilt limit is selected, the full tilt range will be available.		
<input type="checkbox"/> 1223220 Set at 20°	<input type="checkbox"/> 1223230 Set at 30°	<input type="checkbox"/> 1223240 Set at 40°
ⓘ The Stellar GL and GLT offer a Power Tilt option. See separate order forms. pdgmobility.com/power-tilt		

Recline

<input type="checkbox"/> 1223310 Reclining Backrest (Manual)	\$ 895
30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery.	

Front Rigging

<input type="checkbox"/> 1225000 No Front Rigging	<input type="checkbox"/> 1225001 Sector Blocks only	
Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out		
<input type="checkbox"/> 1225410 Fixed 70°		\$ 95
<input type="checkbox"/> 1225420 Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°)	\$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.		
Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out.		
<input type="checkbox"/> 1225130 Swing Away 60°	<input type="checkbox"/> 1225140 Swing Away 70°	
<input type="checkbox"/> 1225200 Swing Away Elevating Legrests	Includes Calf Pads.	E0990 \$ 375
<input type="checkbox"/> 1225300 1" Offset Hangers for Swing Away		\$ 75
<input type="checkbox"/> 1225110 Standard Plug In 70° Front Rigging	Legacy Design. Not compatible with Easy-Out.	

Footplates

<input type="checkbox"/> 1226010 Composite Footplates		
<input type="checkbox"/> 1226020 Aluminum Angle Adjustable Footplates	K0040	\$ 165
<input type="checkbox"/> 1226030 One Piece Fixed Footplate	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/> 1226040 One Piece Flip-up Footplate	Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/> 1226050 Interlocking Angle Adjustable Footplates	Available for Seat Widths up to 21" only.	K0040 \$ 195

Footrest/Legrest Accessories

Calf Pads Not compatible with Gel Padding.		
<input type="checkbox"/> 1227012 For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.	\$ 200
<input type="checkbox"/> 1227010 For Swing Away Front Rigging		\$ 135
<input type="checkbox"/> 1227020 Gel Padding	Not compatible with Calf Pad or Easy-Out Front Rigging.	\$ 125
Gel Knee Pads		
<input type="checkbox"/> 1227032 For Easy-Out™ Angle Adjustable Front Rigging		\$ 200
<input type="checkbox"/> 1227030 For Swing Away ELR Front Rigging		\$ 100
<input type="checkbox"/> 1227040 Adjustable Length Heel Loops		E0951 \$ 50
<input type="checkbox"/> 1227050 Bodypoint Padded Calf Strap	Sizes greater than 22" are not padded.	K0038 \$ 195
Miller's Dynamic Footrest Mounts Not compatible with Easy-Out Front Rigging.		
<input type="checkbox"/> 1227060 For Standard Plug In Front Rigging		\$ 445
<input type="checkbox"/> 1227065 For Swing Away 70° Front Rigging		\$ 575

Front Rigging Extension Tubes


Regular Extension Tubes

<input type="checkbox"/>	1228010	Short	14"-18" for Easy-Out Front Rigging	13"-17" for Swing Away and Standard Plug In Front Rigging		
<input type="checkbox"/>	1228020	Long	18"-24" for Easy-Out Front Rigging	15"-19" for Swing Away and Standard Plug In Front Rigging	\$	75

Extra Strength Extension Tubes (Solid)

Available on Plug In Front Rigging only. Not compatible with Easy-Out Front Rigging.

<input type="checkbox"/>	1228110	Short (13"-17")	<input type="checkbox"/>	1228120	Long (15"-19")	\$	100
--------------------------	---------	-----------------	--------------------------	---------	----------------	----	-----

 **Lower Leg Length** If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.
 Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.

Miscellaneous Options

<input type="checkbox"/>	1230050	Protective Lower Padding			\$	675
<input type="checkbox"/>	1230080	Gel Hip Guides	Reduces available Seat Width by 2".	E0956 & E1028	\$	395
<input type="checkbox"/>	1230110	IV Pole	Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	E1029	\$	295
<input type="checkbox"/>	1230120	O2 Holder	Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	E2208	\$	395
<input type="checkbox"/>	1230170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.		\$	10

On-Chair Seating & Positioning

 Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions.
 Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.